

# Seclira®

## Gel Cockroach Bait

# The gel partner you have been waiting for

**Seclira® Gel Cockroach Bait** has been specifically designed as a crevice, crack or spot treatment for the control of cockroaches. This highly palatable formulation is powered by a fast-acting active ingredient and BASF Genuine Transfer Effect™ technology providing rapid, whole of colony control. Seclira Gel Cockroach Bait also combines with the suite of Seclira products and partners in rotation with Goliath® Gel Cockroach Bait for unmatched performance in all conditions.

- Highly palatable formulation backed by BASF Genuine Transfer Effect Technology for entire colony control
- Unique, fast acting active ingredient for rapid control
- Developed to combine with the full Seclira range of products
- Suitable for commercial and domestic use, Seclira Cockroach Gel Bait is your perfect rotation partner with Goliath Cockroach Gel

For more information on  
Seclira® Gel Cockroach Bait  
visit [pest-control.basf.com.au](https://pest-control.basf.com.au)  
or call **1800 558 399**

**BASF**  
We create chemistry

# Seclira®

## Gel Cockroach Bait

Seclira Gel Cockroach Bait is a specifically designed spot treatment for the targeted control of German and other common cockroaches.

Its strong palatability, appealing to gel-bait averse cockroaches, provides a solution to resistance management with a strong transfer effect. Entire colony control is achieved with minimal behaviour, eating and social changes, providing the foundation for its strong, long lasting residual. Available in a convenient 30g pack size, Seclira Gel Cockroach Bait is suitable for domestic and commercial use, with a targeted treatment in areas of known cockroach activity.

SITUATION	PEST	RATES	REMARKS
Domestic, Commercial and Public Service Buildings	German and other common cockroaches	Apply 1-3 spots per m <sup>2</sup> . Recommended spot size is 0.25-0.5 g. Increase the number of spots (up to 3) and spot size depending upon severity of infestation and species present.	Gel placement should be at or near areas where cockroaches gather, such as in cracks and crevices

Refer to the product label for full directions and use table.



## Packaging

Packaging – 4 x 30g Syringes

For more information on Seclira Gel Cockroach Bait visit [pest-control.basf.com.au](http://pest-control.basf.com.au) or contact your local BASF representative on **1800 558 399**

### ALWAYS READ AND FOLLOW LABEL DIRECTIONS BEFORE USING ANY PRODUCT IN THIS FACT SHEET.

This fact sheet is intended as general advice. Disclaimer: The information submitted in this publication is based on current BASF knowledge and experience. In view of the many factors that may affect its application, this data does not relieve the user from carrying out their own tests. The data does not imply assurance of certain properties or of suitability for a specific purpose. It is the responsibility of the user to ensure that any proprietary rights and existing laws and legislation are observed.